

# Greencare<sup>®</sup>

MEASURED MOBILITY



# easy1<sup>®</sup>



A lightweight family of chairs with a carry weight with wheels and footrests removed at around 7kg that offer self-propelling or carer controlled operation. With a user weight of 125kg these chairs are both strong and yet remarkably lightweight. They are also safe as they have been tested to British and International Standards ISO 7176-19 which includes successfully passing crash testing making them safe to use in transport.

These comfortable chairs offer the option of brakes on both self-propelling and carer-controlled options with quick release wheels on the self-propelling chair as standard. There is an option of a tension adjustable back to further improve comfort.

## FEATURES

- Seat widths 16", 17", 18" and 20" available for comfort and super strength.
- Seat depths 17".
- Backrest and seat are made in cushioned fabric for greater comfort with tension adjustment on back canvas for additional contoured support.
- A lap belt is fitted as standard
- A calf strap is provided
- Greencare wheelchairs are crash tested to ISO 7176-19.



YOUR LOCAL DEALER



SWING AWAY  
ARMRESTS



TENSION  
ADJUSTABLE BACK



ELEVATING  
LEG REST



BRAKES TO ASSIST  
CARERS PUSHING  
THE CHAIR



FOLDING  
BACK POSTS

To further reduce chair size



ALL CHAIRS FOLD

Making it easier for  
transportation and storage.

### CUSHIONS



Viscoe elastic and high resilience foam core washable with non slip on base.